

---

# THORLABS

## MDFM-MF2

# Olympus AX, BX, and IX Microscopy Filter Cube

Assembly Manual



---

# Table of Contents

---

<i>Part 1. Safety</i> .....	<i>3</i>
<i>Part 2. Assembly Instructions</i> .....	<i>3</i>
<i>Part 3. Regulatory</i> .....	<i>9</i>
3.1. Waste Treatment is Your Own Responsibility .....	<i>9</i>
3.2. Ecological Background .....	<i>9</i>
<i>Part 4. Thorlabs' Worldwide Contacts</i> .....	<i>10</i>

## Part 1. Safety

All statements regarding safety of operation and technical data in this instruction manual will only apply when the unit is operated correctly. Handle all optics with care. Gloves are recommended for optic installations. Take care when tightening optics, as overtightening optics may cause the optic to break or shatter.

## Part 2. Assembly Instructions

### Step 1

Carefully unpack the Olympus MDFM-MF2 Filter Cube assembly. Take note where the emitter and exciter sides of the cube are (see Figure 1).




Figure 1. MDFM-MF2 Filter Cube Assembly

Note the orientation of the retainer rings. The flat side of both retaining rings should always face away from the cube. The side of the retaining ring with the mounting lip should face inwards towards the filter and dichroic mirror.

### Step 2

Remove the two screws holding the filter cube together. Carefully separate the two filter cube halves.





Figure 2. Disassembling the filter cube

### Step 3

Remove the emitter side of the filter cube and orient so the emitter port is facing down as shown in Figure 3. Loosen the two screws holding the dichroic mirror retainer clip. Do not remove the screws.



**Figure 3. Loosening the dichroic mirror retainer clip**


### Step 4

The MDFM-MF2 Filter Cube holds a standard set of Thorlabs MDF Series Dichroic Mirror, Emission, and Excitation Filter Set, sold separately. The MDF-CY3.5 Filter set is shown in Figure 4.



**Figure 4. Dichroic mirror set**

The emission and excitation filters can be easily differentiated by the thickness of the filter mount rings: emission (3.5 mm) and excitation (5.0 mm). Note that the dichroic mirror has a coated and an uncoated side (see Figure 5).



**Figure 5. Dichroic mirror orientation**

When held at an angle, the coated side will appear uniform in reflection all the way to the edge. The uncoated side, when held at the same angle, will not have a uniform reflection near the edge. There will be a band of transmission that appears through the uncoated substrate. It is important to determine the coated and uncoated sides of the dichroic mirror as there is only one orientation that will work with the Filter Cube.

## Step 5

Place the dichroic mirror under the retaining clip, making sure the mirror seats properly in the recess. The reflective surface of the dichroic mirror should face **UP**.



**Figure 6. The dichroic mirror mounted in the emitter half cube**

Slowly tighten the retaining clip screws until snug. Caution: Do not over tighten the retaining clip screws as this may damage the dichroic mirror.

## Step 8

Reassemble the filter cube by placing the half of the assembly with the mounted dichroic mirror into the empty half. Align the two halves using the two assembly screw holes. Slowly tighten the assembly screws, alternating between them, until snug. Caution: Do not over tighten the assembly screws, as doing so may strip the threads or screw heads. At this point the dichroic mirror should be facing the exciter side of the filter cube (see Figure 7).




Figure 7. Looking into the exciter side of the filter cube

## Step 9

Orient the filter cube so the emitter side is facing up. Place the emission filter into the mount with the filter arrow facing up, away from the mount. Carefully thread the retainer ring (lip side down) onto the cube.





Figure 8. Installing the emission filter



Note: The threads are not parallel with the face of the cube. Take extra care to make sure the retainer ring is threaded properly onto the cube. The ring should rotate smoothly and easily when threaded properly. Tighten the retainer ring using the supplied spanner wrench, or Thorlabs Spanner, until snug.




**Figure 9.** Note the orientation of the emission filter and retainer ring

## Step 10

Place the excitation filter into the filter cube with the arrow facing down into the mount. Carefully thread the retainer ring (lip side facing inwards) onto the cube.

Take extra care to make sure the retainer ring is threaded properly onto the cube. The ring should rotate smoothly and easily when threaded properly.

Tighten the retainer ring using the supplied spanner wrench.



**Figure 10.** Installing the excitation filter

**Step 11**

The completed MDFM-MF2 Olympus Filter Cube is shown in Figure 11.



*Figure 11. The MDFM-MF2 Filter Cube with emission filter, excitation filter, and dichroic mirror installed*



## Part 3. Regulatory

As required by the WEEE (Waste Electrical and Electronic Equipment Directive) of the European Community and the corresponding national laws, Thorlabs offers all end users in the EC the possibility to return “end of life” units without incurring disposal charges.

- This offer is valid for Thorlabs electrical and electronic equipment:
- Sold after August 13, 2005
- Marked correspondingly with the crossed out “wheelie bin” logo (see right)
- Sold to a company or institute within the EC
- Currently owned by a company or institute within the EC



*Wheelie Bin Logo*

- Still complete, not disassembled and not contaminated

As the WEEE directive applies to self-contained operational electrical and electronic products, this end of life take back service does not refer to other Thorlabs products, such as:

- Pure OEM products, that means assemblies to be built into a unit by the user (e. g. OEM laser driver cards)
- Components
- Mechanics and optics
- Left over parts of units disassembled by the user (PCB's, housings etc.).

If you wish to return a Thorlabs unit for waste recovery, please contact Thorlabs or your nearest dealer for further information.

### **3.1. Waste Treatment is Your Own Responsibility**

If you do not return an “end of life” unit to Thorlabs, you must hand it to a company specialized in waste recovery. Do not dispose of the unit in a litter bin or at a public waste disposal site.

### **3.2. Ecological Background**

It is well known that WEEE pollutes the environment by releasing toxic products during decomposition. The aim of the European RoHS directive is to reduce the content of toxic substances in electronic products in the future.

The intent of the WEEE directive is to enforce the recycling of WEEE. A controlled recycling of end of life products will thereby avoid negative impacts on the environment.

## Part 4. Thorlabs Worldwide Contacts

### USA, Canada, and South America

Thorlabs, Inc.  
435 Route 206  
Newton, NJ 07860  
USA  
Tel: 973-579-7227  
Fax: 973-300-3600  
www.thorlabs.com  
email: feedback@thorlabs.com

### Europe

Thorlabs GmbH  
Hans-Böckler-Str. 6  
85221 Dachau  
Germany  
Tel: +49-(0)8131-5956-0  
Fax: +49-(0)8131-5956-99  
www.thorlabs.de  
email: Europe@thorlabs.com

### France

Thorlabs SAS  
109, rue des Côtes  
78600 Maisons-Laffitte  
France  
Tel: +33 (0) 970 444 844  
Fax: +33 (0) 811 381 748  
www.thorlabs.de  
email: slaes.fr@thorlabs.com

### Japan

Thorlabs Japan Inc.  
Higashi Ikebukuro  
Q Building 1<sup>st</sup> Floor 2-23-2  
Toshima-ku, Tokyo 170-0013  
Japan  
Tel: +81-3-5979-8889  
Fax: +81-3-5979-7285  
www.thorlabs.jp  
email: sales@thorlabs.jp

### UK and Ireland

Thorlabs LTD.  
1 Saint Thomas Place, Ely  
Cambridgeshire CB7 4EX  
Great Britain  
Tel: +44 (0)1353-654440  
Fax: +44 (0)1353-654444  
www.thorlabs.de  
email: sales.uk@thorlabs.com

### Scandinavia

Thorlabs Sweden AB  
Box 141 94  
400 20 Göteborg  
Sweden  
Tel: +46-31-733-30-00  
Fax: +46-31-703-40-45  
www.thorlabs.de  
email: scandinavia@thorlabs.com

### China

Thorlabs China  
Oasis Middlering Centre  
3 Building 712 Room  
915 Zhen Bei Road  
Shanghai  
China  
Tel: +86 (0)21-32513486  
Fax: +86 (0)21-32513480  
www.thorlabs.com  
email: chinasales@thorlabs.com



Thorlabs, Inc.  
435 Route 206N  
Newton, NJ 07860 USA

Phone: (973) 579-7227 ♦ Fax: (973) 300-3600  
[www.thorlabs.com](http://www.thorlabs.com)

---